

AccuStart™ MarkerSelect PCR SuperMix ROX™

Cat. No.	95104-250	Size:	250 x 20- μ L reactions (2 x 1.25 mL)
	95104-012		1250 x 20- μ L reactions (10 x 1.25 mL)
	95104-05K		5000 x 20- μ L reactions (1 x 50 mL)

Store at **-20°C protected from light**

Description

AccuStart MarkerSelect PCR SuperMix is a 2X concentrated, ready-to-use reaction cocktail for routine PCR amplification of DNA fragments for characterization with allelic discrimination assays. It contains all components, except primers, probes and template. AccuStart MarkerSelect PCR SuperMix simplifies reaction assembly, improves reproducibility, and reduces the risk of contamination. A key component is AccuStart™ Taq DNA polymerase, which contains monoclonal antibodies that bind to the polymerase and keep it inactive prior to the initial PCR denaturation step. Upon heat activation the antibodies denature irreversibly enabling specific and efficient primer extension and probe utilization with the convenience of room temperature reaction assembly. AccuStart MarkerSelect mixes also contain additives to enable robust and accurate genotyping from agricultural and environmental samples that might contain PCR inhibitors.

Components

AccuStart MarkerSelect PCR SuperMix 2X reaction buffer containing MgCl₂, dNTPs (dATP, dCTP, dGTP, dTTP), AccuStart Taq DNA Polymerase, ROX reference dye, and stabilizers.

Storage and Stability

AccuStart MarkerSelect PCR SuperMix is stable for 1 year when stored in a constant temperature freezer at -20°C and protected from light. For convenience, it may be stored unfrozen at +2 to +8°C for up to 6 months. No loss of performance was demonstrated after 20 freeze-thaw cycles or 2 months at +20°C.

Reaction Assembly

Component	Volume for 20- μ L rxn.	Final Concentration
AccuStart MarkerSelect PCR SuperMix ROX (2X)	10 μ L	1X
Primer & Probe	Variable	0.5X to 1X
Nuclease-free water	Variable	
DNA Template	0.5 – 5 μ L	Variable
Final Volume (μ L)	20 μ L	

Note: For smaller reaction volumes (i.e. 5- μ L reactions), scale all components proportionally.

Reaction Protocol

Incubate the completed reaction mix in thermal cycler as follows:

Initial denaturation:	95°C, 5 min (up to 10 min for difficult templates)
PCR cycling (20 – 40 cycles):	95°C, 5 s 60°C, 30s
Hold	4°C until processed for analysis

Full activation of AccuStart Taq DNA polymerase occurs within 30 seconds at 95°C. Complete denaturation of dsDNA target is important for efficient PCR amplification and may require different initial denaturation times depending on the properties of a given target sequence. Cycling protocol may require adjustment for longer amplicons (>150bp).

Quality Control

Kit components are free of contaminating DNase and RNase and validated in PCR performance tests.

Limited Label Licenses

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,804,375, 5,538,848, 5,723,591, 5,876,930, 6,030,787 and 6,258,569. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained from the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

The purchase of this product includes a limited, non-transferable right to use the purchased amount of the product to perform Applied Biosystems' patented Passive Reference Method for the purchaser's own internal research. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. For information about these rights or on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

Licensed to Quanta BioSciences, under U.S. Patent Nos. 5,338,671, 5,587,287, and foreign equivalents for use in research only..

AccuStart is a trademark of Quanta BioSciences Inc. ROX is a trademark of Life Technologies Corporation.

©2010 Quanta BioSciences, Inc. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.