



PerfeCta® Multiplex qPCR SuperMix, Low ROX™

Cat. No. 95108-050 Size: 50 x 50- μ L reactions
95108-200 200 x 50- μ L reactions
95108-01K 1000 x 50- μ L reactions

Store at -20°C protected from light

Description

PerfeCta Multiplex qPCR SuperMix, Low ROX is a 2X concentrated, ready-to-use reaction cocktail that contains all components, except primers, probe(s), and template for real-time quantitative PCR on Applied Biosystems 7500, 7500 Fast, ViiA7 or Stratagene MX series of real-time PCR systems. The system transcends multiplex limitations of conventional PCR master mixes, enabling unbiased amplification of up to five target sequences in a single tube. Suppression of low copy amplicons by high copy reference targets in the amplification is a common problem in multiplex PCR. This can skew, or mask the apparent representation and quantification of low copy target sequences. PerfeCta Multiplex qPCR SuperMix, Low ROX delivers dynamic range and sensitivity to multiplexed qPCR that is comparable to that for singleplex qPCR probe assays without the need for limiting or variable primer concentrations.

The robust qPCR performance of PerfeCta Multiplex qPCR SuperMix, Low ROX also overcomes suppression by common PCR inhibitors in blood or environmental samples, and is a highly sensitive reagent for quantification of low copy target sequences in challenging samples.

A key component of this supermix is AccuStart™ Taq DNA polymerase, which contains monoclonal antibodies that bind to the polymerase and keep it inactive prior to the initial PCR denaturation step. Upon heat activation (2 minutes at 95°C), the antibodies denature irreversibly, releasing fully active, unmodified Taq DNA polymerase. This enables specific and efficient primer extension with the convenience of room temperature reaction assembly.

Instrument Compatibility

Different real-time PCR systems employ different strategies for the normalization of fluorescent signals and correction of well-to-well optical variations. It is critical to match the appropriate qPCR reagent to your specific instrument. PerfeCta Multiplex qPCR SuperMix, Low ROX provides seamless integration on the Applied Biosystems 7500, 7500 Fast, ViiA™7 or Stratagene MX series of real-time PCR systems. These ©2016 Quanta BioSciences, Inc. All rights reserved.

Intended for molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease.

95108 / IFU-074.1 REV 1

Instruments utilize variable excitation wavelengths that are tuned to the each respective dye detection channel and provide superior sensitivities and dynamic ranges for multiplex probe applications. Your choice of probe reporter dyes and any optional internal reference dye must be matched to the excitation and emission optics of your particular instrument. Please consult the user manual for your real-time PCR system.

Components

PerfeCta Multiplex qPCR Supermix, Low ROX 2X reaction buffer containing optimized concentrations of MgCl₂, dNTPs (dATP, dCTP, dGTP, dTTP), AccuStart Taq DNA Polymerase, ROX Reference Dye (for 580-585 nm excitation), and stabilizers.

Storage and Stability

PerfeCta Multiplex qPCR SuperMix, Low ROX is stable for 1 year when stored in a constant temperature freezer at -20°C, protected from light. For convenience, it may be stored unfrozen at +2 to +8°C for up to 6 months. Repeated freezing and thawing of the supermix does not affect performance.

Guidelines for Multiplex qPCR:

- The design of highly specific primers and probes is a critical and challenging aspect of successful multiplex qPCR. Each primer and probe should have similar thermodynamic properties to support efficient PCR amplification using a common temperature cycling program for all amplicons. The use of computer aided primer design programs is encouraged in order to minimize the potential for internal secondary structure and complementation at 3'-ends within each primer, primer pairs, and primer/probe combinations
- Amplicon size should be consistent for each target sequence and limited to approximately 65 - 100 bp.
- Limiting primer concentration for high copy genes is acceptable, but not required. A final concentration of 300 nM each primer and 100 to 250 nM probe is effective for most applications. Each probe for a multiplex assay should be labeled using dyes with minimal spectral overlap and non-fluorescent quencher compounds. Matching dyes with discrete fluorescent excitation and emission optima improves the accuracy of the multicomponenting, or dye deconvolution algorithms employed by the real-time PCR analysis software.
- Preparation of a reaction cocktail is recommended to reduce pipetting errors and maximize assay precision. Assemble the reaction cocktail with all required components except sample template (genomic DNA or cDNA) and dispense equal aliquots into each reaction tube. Add the DNA template to each reaction as the final step. Addition of samples as 5 to 10- μ L volumes will improve assay precision.

©2016 Quanta BioSciences, Inc. All rights reserved.

Intended for molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease.

95108 / IFU-074.1 REV A

Guidelines for qPCR continued:

- Suggested input quantities of template are: cDNA corresponding to 10 pg to 1 µg of total RNA; 100 pg to 1 µg genomic DNA
- After sealing each reaction, vortex gently to mix contents. Centrifuge briefly to collect components at the bottom of the reaction tube.

Reaction Assembly

| Component | Volume for 50-µL rxn. | Final Concentration |
|--|-----------------------|---------------------|
| PerfeCt _a Multiplex qPCR SuperMix | 25 µL | 1x |
| Forward primers | variable | 100 – 500 nM each |
| Reverse primers | variable | 100 – 500 nM each |
| Probes | variable | 100 – 250 nM each |
| Nuclease-free water | variable | |
| Template(s) | <u>5 – 10 µL</u> | variable |
| Final Volume (µL) | 50 µL | |

Note: Reaction volume can be scaled from 10 to 50 µL depending on the reaction plate (i.e. 384-well vs. 96-well) and qPCR system. Scale all components proportionally to the desired final reaction volume. When performing qPCR in 384-well plates, use reaction volumes between 10 and 15 µL. High (20-µL) volumes can impair mixing of sample template with the reaction cocktail and can compromise qPCR performance and reproducibility

Cycling Protocol

Incubate complete reaction mix in a real-time thermal detection system as follows:

| | Standard Protocol | Fast Protocol |
|---------------------------------|---|--|
| Initial denaturation/activation | 95°C, 2 to 3 min | 95°C, 30s |
| PCR Cycling (35 to 45 cycles): | 95°C, 10 to 15 s 55 – 65°C, 60s (collect and analyze data) | 95°C, 5s 60°C, 30s (collect and analyze data) |

Full activation of AccuStart Taq DNA polymerase occurs within 30 seconds at 95°C; however, optimal initial denaturation time is template dependent and will affect qPCR efficiency and sensitivity. Amplification of genomic DNA targets benefit from a prolonged initial denaturation step (5-10 min) to fully denature and fragment the template. This minimizes the potential for renaturation of long fragments and/or repetitive sequence regions that can impair replication of the target sequence by the PCR process.

Quality Control

Kit components are free of contaminating DNase and RNase. PerfeCt_a Multiplex qPCR SuperMix, Low ROX is functionally tested in a four-plex TaqMan qPCR using variable concentrations of one target sequence from 100 to 1 x 10⁷ copies and 1 x 10⁸ copies each of three other target sequences. Kinetic analysis must demonstrate linear resolution over six orders of dynamic range ($r^2 > 0.995$) and a PCR efficiency > 90%.

Limited Label Licenses

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,804,375, 5,538,848, 5,723,591, 5,876,930, 6,030,787 and 6,258,569. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained from the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

The purchase of this product includes a limited, non-transferable right to use the purchased amount of the product to perform Applied Biosystems' patented Passive Reference Method for the purchaser's own internal research. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. For information about these rights or on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

Licensed to Quanta BioSciences, under U.S. Patent Nos. 5,338,671, 5,587,287, and foreign equivalents for use in research only.

PerfeCt_a and AccuStart are trademarks of Quanta BioSciences Inc. TaqMan is a registered trademark of Roche Molecular Systems, Inc. Applied Biosystems, ViiA and ROX are trademarks of Life Technologies Corporation. Stratagene, MX3000P, MX3005P and MX4000 are trademarks of Agilent Technologies, Inc.