



10mM dNTP Mix

Cat. No 95062-200 Size: 200 μ L (2 pmol each dATP, dCTP, dGTP, dTTP)
95062-01K 1000 μ L (10 pmol each dATP, dCTP, dGTP, dTTP)

Store at **-25°C to -15°C**

Description

The 10 mM dNTP Mix is a solution of high purity deoxyribonucleoside triphosphates that has been functionally qualified for real-time quantitative PCR (qPCR). It is suitable for use in conventional end-point PCR, real-time qPCR, first-strand cDNA synthesis, as well as other applications that require dNTP as substrate.

We have found that dNTP quality and purity can directly affect dynamic range, sensitivity and precision of qPCR assays. Surprisingly, dNTP preparations that fail our qPCR functional testing are often >99% triphosphate purity by HPLC analysis and function perfectly well in conventional end-point PCR assays, including amplification of long (up to 20 kb) templates. The use of qPCR certified dNTPs is a simple but important measure to ensure consistent assay performance.

Storage and Stability

Store components in a constant temperature freezer at -25°C to -15°C upon receipt.

Repeated freezing and thawing of the reaction mix is not recommended.

For lot specific expiry date, refer to package label, Certificate of Analysis or Product Specification Form.

Components

10 mM dATP, 10 mM dCTP, 10 mM dGTP, and 10 mM dTTP

Quality control

10 mM dNTP mix is functionally tested in real-time PCR for amplification of a single copy gene from human genomic DNA and is free of detectable RNase or DNase. Each deoxyribonucleoside triphosphate is > 99% pure by HPLC analysis.

Limited Label Licenses

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this manual and for use with components contained in the kit only. QIAGEN Beverly, Inc. grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the protocols provided with the product, this manual, and additional protocols available at www.quantabio.com. Some of these additional protocols have been provided by Quantabio product users. These protocols have not been thoroughly tested or optimized by QIAGEN Beverly, Inc.. QIAGEN Beverly, Inc. neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN Beverly, Inc. makes no warranty that this kit and/or its use(s) do not infringe the rights of third-parties.
3. This kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN Beverly, Inc. specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN Beverly, Inc. may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

©2018 QIAGEN Beverly Inc. 100 Cummings Center Suite 407J Beverly, MA 01915

Quantabio brand products are manufactured by QIAGEN, Beverly Inc.

Intended for molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease.