



## Extracta™ DBS

Cat No: 95171-010  
95171-500

Size: 10 mL  
500 mL

Store at 4 °C

### Description

Extracta DBS is a ready-to-use DNA extraction reagent for rapid and efficient recovery of PCR-ready DNA from dried blood spots (DBS) on Guthrie cards or Whatman 903 filter paper. This patented single-solution process produces DNA that is suitable for a variety of downstream applications including real-time qPCR or template generation for Next Generation Sequencing (NGS) or Sanger sequencing (1-3). DNA eluates are substantially free of PCR inhibitors and compatible with a variety of end-point PCR or real-time qPCR reagents. Application of Extracta DBS with PerfeCTa® qPCR ToughMix® or PerfeCTa® MultiPlex qPCR ToughMix enables accurate and reproducible quantification of DNA sequences in blood using TaqMan® hydrolysis probe real-time qPCR.

### Components

Extracta DBS 1X concentrated, ready-to-use, proprietary DNA extraction reagent

### Storage and Stability

Extracta DBS is stable for up to 1 year at 2-8°C. Please consult the product box label for expiration date information.

### Recommended Protocol for DNA Isolation from Dried Blood Spots (DBS)

1. Place a single 3.2 mm punch from DBS specimen in each PCR tube or well of a PCR plate.
2. Add 100 µL of Extracta DBS reagent and rinse the DBS punch by gentle vortexing or pipetting up and down.
3. Cap the tubes (or seal the plate) and centrifuge at ~3,500 rpm for 5 minutes.
4. Using a pipettor, remove as much liquid and debris as possible and discard the rinse into 10% bleach or in accordance with your institution's safety policies for processing dried blood spot specimens.
5. Add 50 µL of fresh Extracta DBS reagent to each rinsed DBS specimen.
6. Cap the tubes (or seal the plate) and then briefly centrifuge (10-30s) at ~3,500 rpm to ensure that the punch is completely submerged in the reagent.
7. Place the tubes or plate in a thermal cycler with a heated lid and incubate at 95°C for 20 minutes, then hold at 4-8°C.
8. Remove the tubes or plate from the thermal cycler and use up to 5 µL of extracted DNA as template in each 20-25-µL PCR.

Note: it is normal for the extracted DNA solution to have a pink color.

Extracted DNA should be used immediately in down-stream amplification procedures. It may be refrigerated at 2-8°C for short term storage before PCR amplification. The effect of prolonged storage on specific applications may vary.

### Reagents and Equipment Required but not Provided

- Microcentrifuge tubes (0.5 mL or 1.5 mL), PCR tubes (0.2 mL) or multiwell plates
- P100 / P200 pipettor or multichannel pipettor
- PCR reagents
- Thermal cycler

©2017 Quanta BioSciences, Inc. All rights reserved.

Intended for molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease.

95171 / IFU-113.1 REVA

## Related Products

PerfeCTa qPCR ToughMix, 250R Cat. No. 95112-250,  
PerfeCTa qPCR ToughMix, 5000R, Cat. No. 95112-05K,  
PerfeCTa qPCR ToughMix, 1250R, Cat. No. 95112-012

PerfeCTa qPCR ToughMix ROX, 250R, Cat. No. 95113-250,  
PerfeCTa qPCR ToughMix ROX, 5000R, Cat. No.95113-05K,  
PerfeCTa qPCR ToughMix ROX, 1250R, Cat. No. 95113-012

PerfeCTa qPCR ToughMix L-ROX, 250R, Cat. No.95114-250,  
PerfeCTa qPCR ToughMix L-ROX, 5000R, Cat. No.95114-05K,  
PerfeCTa qPCR ToughMix L-ROX, 1250R, Cat. No.95114-012

PerfeCTa Multiplex qPCR ToughMix, 250R, Cat. No.95147-250,  
PerfeCTa Multiplex qPCR ToughMix, 1000R, Cat. No.95147-01K,  
PerfeCTa Multiplex qPCR ToughMix, 5000R, Cat. No.95147-05K

PerfeCTa Multiplex qPCR ToughMix, ROX 250R, Cat. No. 95148-250,  
PerfeCTa Multiplex qPCR ToughMix, ROX 1000R, Cat. No. 95148-01K,  
PerfeCTa Multiplex qPCR ToughMix, ROX 5000R, Cat. No. 95148-05K

PerfeCTa Multiplex qPCR ToughMix, L-ROX 250R, Cat. No. 95149-250,  
PerfeCTa Multiplex qPCR ToughMix, L-ROX 1000R, Cat. No.95149-01K,  
PerfeCTa Multiplex qPCR ToughMix, L-ROX 5000R, Cat. No.95149-05K

## Limited Label License

This product is licensed under US Patent 9,206,468 and foreign equivalents for non-commercial research purposes only.

## References

- 1) Baker et al. 2016. Improving newborn screening for cystic fibrosis using next-generation sequencing technology: a technical feasibility study. Genet Med. 18(3):231-238.
- 2) Baker et al. 2010. Implementing Routine Testing for Severe Combined Immunodeficiency within Wisconsin's Newborn Screening Program. Public Health Rep. 125 Suppl 2, 88-95.
- 3) Verbsky et al. 2012. Newborn Screening for Severe Combined Immunodeficiency; The Wisconsin Experience (2008-2011). J. Clin. Immunol. 32, 82-88.

## Trouble Shooting Guide

| Problem   | Possible Cause   | Solution   |
|---|--|--|
| No PCR product or non-specific product (or signal) from positive control samples using purified genomic DNA | PCR primers, reagents or cycling conditions were not optimal | Refer to the appropriate PCR reagent product manual to optimize PCR conditions   |
| No PCR product or non-specific product (or signal) from extracts  | Inadequate extract heating                                   | Ensure that DBS punch extracts are incubated at $\geq 95^{\circ}\text{C}$ . Program thermal cycler for $96^{\circ}\text{C}$ .  |
|   | Extraction time was too short                                | Incubate sample in extraction reagent for up to 30 minutes at $95^{\circ}\text{C}$ .   |
|   | Too much extract or template in PCR                          | Use less than 1/10 volume of extract in the PCR reaction. Extracts can be diluted 5-, 10-, 20-fold or more in TE buffer (10 mM Tris-HCl pH 8.0, 0.1 mM EDTA) prior to PCR. |

Extracta, PerfeCTa, and ToughMix are trademarks of Quanta BioSciences Inc. TaqMan is a registered trademark of Roche Molecular Systems, Inc.

©2017 Quanta BioSciences, Inc. All rights reserved.

Intended for molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease.

95171 / IFU-113.1 REVA