


Rapid and accurate quantification of Illumina NGS libraries using the Q real-time qPCR Instrument

Achieve 60% faster time to reliable results



Peter Bartholomew, Hongbo Liu, Eleanor Kolossovski, Brian Komorous, and David Schuster
100 Cummings Centers 407J, Beverly, MA. 01915

Introduction

Accurate quantification of the number of amplifiable library molecules is a critical factor for obtaining high quality read data with next-generation sequencing technologies. The high sensitivity, broad dynamic range, and specificity of qPCR to quantify library molecules that are suitable for the bridge PCR provide significant advantages over methods for total DNA quantification. However, these advantages are often offset by the time to result, requirement for inclusion of absolute DNA standards in every qPCR run, and errors associated with dilution of libraries so that reportable results are within the linear dynamic range of the technology. Here we describe application of a new real-time quantitative PCR instrument, the Q from QuantaBio, to simplify reliable library quantification with faster run times.








Features of the Q Real-Time Quantitative PCR Instrument

- Magnetic Induction Technology**
Rapidly heats reactions held in a unique spinning aluminum rotor
- Superior Temperature Uniformity of ± 0.05°C**
Eliminates well position effects associated with traditional peltier block-based real time cyclers
- Ultra-Fast Data Acquisition**
Robust, fixed optical path allows simultaneous acquisition of all channels with no need for reference dyes or crosstalk compensation
- Scalable and Wireless**
Up to 10 Q instruments can be operated from a single workstation wirelessly via Bluetooth, enabling processing of 480 samples simultaneously
- Portable and Compact**
The compact size and 4.5 pound weight of the Q allows easy portability with no need for calibration
Compact size and 4.5 pound weight of the Q allows easy portability with no need for calibration while occupying 1/4 the bench space of traditional cyclers
- Powerful Software**
User-friendly Q-qPCR software for advanced automated statistical analysis including relative quantification, absolute quantification, genotyping and allelic discrimination

Materials and Methods

Most trials were conducted on a single Q using DNA standards, primers and SYBR Green Supermix from the PerfeCTa NGS Library Quantification Kit for Illumina Sequencing Platforms (QuantaBio cat# 95154). PerfeCTa SYBR Green Fastmix (QuantaBio cat# 95072) was also used where noted. The 426-bp unknown test sample was prepared by pooling an arbitrary amount of each QuantaBio DNA standard. DNA libraries were prepared from different microbial DNA sources using the sparQ DNA Frag & Library Prep Kit (QuantaBio cat# 95194). Multiple library samples representing a range of GC-contents were pooled and analyzed by Agilent Bioanalyzer to establish the average fragment size value of 450-bp. Data was analyzed using the Q-qPCR software (v1.0.0).

Results



Conclusions

Together, the results presented clearly establish the suitability of the Q real-time PCR instrument and PerfeCTa NGS Quantification Kit for quantification of NGS libraries of various sizes and GC contents.

The clear benefits provided by the Q for NGS Quantification include:

- Highly precise measurements across multiple trials
- High efficiency amplifications under varied cycling conditions
- Exceptional quantitative sensitivity for distinguishing down to 1.2-fold differences
- Reliable results and performance from run times 60% shorter than typical cycling protocols