

## repliQa HiFi ToughMix®

Cat No.	95200-025	Size:	25 x 50 µL reactions (1 x 0.625 mL)
	95200-100		100 x 50 µL reactions (1 x 2.50 mL)
	95200-500		500 x 50 µL reactions (1 x 12.50 mL)

Store at -25°C to -15°C

### Description

The repliQa HiFi ToughMix is a unique, next generation 2x master mix that has 90x higher fidelity compared to *Taq* polymerase. The ToughMix has extreme speed, with extension times as fast as 1-10 sec/kb depending on target length. Additionally, the ToughMix has long range amplification properties as it can amplify fragments up to 24 kb from complex genomic DNA templates or fragments up to 40 kb from virus DNA templates such as *Escherichia virus Lambda* DNA.

The ToughMix is formulated with a genetically modified DNA polymerase coupled with hot start antibodies. It has 5' → 3' polymerase activity, 3' → 5' exonuclease activity, and generates blunt-ended products while providing the ability to amplify through uracils and primers containing inosine or uracil. It is *Tough Tested*, and is tolerant to multiple PCR inhibitors.

### Components

<b>repliQa HiFi ToughMix</b>	2x reaction buffer containing optimized concentrations of MgCl <sub>2</sub> , dNTPs and proprietary formulated HiFi polymerase, hot start antibodies and ToughMix chemistry
------------------------------	---

### Storage and Stability

Store kit components in a constant temperature freezer at -25°C to -15°C protected from light upon receipt. For lot specific expiry date, refer to package label, Certificate of Analysis or Product Specification Form.

### Guidelines for PCR

- The design of highly specific primers is a critical parameter for successful PCR. The use of computer aided primer design programs is encouraged in order to minimize the potential for internal secondary structure. For best results, primer size should be limited to 22 - 35 bp with a melting point of at least 63°C. Ideal GC-content of the primers is 45-60%. A final concentration of 300 nM each primer is effective for most applications. Primers with inosine (dI) and uracil (dU) are acceptable.
- Preparation of a reaction cocktail is recommended to reduce pipetting errors and maximize assay precision. Assemble the reaction cocktail with all required components except sample template and dispense equal aliquots into each reaction tube. Add the DNA template to each reaction as the final step. Addition of samples as 2 to 5-µL volumes will improve assay precision.
- Suggested input quantities of template are: genomic DNA ≤ 200 ng; plasmid DNA ≤ 50 ng; cDNA ≤ 750 ng.
- After sealing each reaction, vortex gently to mix contents. Centrifuge briefly to collect components at the bottom of the reaction tube.
- Longer targets may require a higher primer melting temperature of at least 65°C, and a lower primer concentration of 150 nM.

### Reaction Assembly

Component	Volume for 50 µL rxn	Final Concentration
repliQa HiFi ToughMix (2X)	25 µL	1x
Forward primer	variable	300 nM
Reverse primer	variable	300 nM
Nuclease-free water	variable	
Template	2 – 5 µL	variable
Final Volume (µL)	50 µL	

**Note:** For smaller or larger reaction volumes scale all components proportionally.



## PCR Cycling Protocol

Initial denaturation and final extension steps are typically not required with repliQa HiFi ToughMix.

The 10 s at 98°C during cycling are sufficient to fully activate the HotStart mechanism. For longer fragments, (>10 kb), an initial denaturation of 98°C for 30 s can be added to facilitate denaturation of the DNA template.

PCR cycling (25 - 45 cycles):

\*We recommend using 2-step cycling first. If you are having difficulty, we suggest trying 3-step cycling.

2-Step Cycling*	3-Step Cycling
98°C, 10 s	98°C, 10 s
	(T <sub>m</sub> -5) °C, 5 s
68°C	68°C
≤ 1 kb: 1 sec	≤ 1 kb: 1 sec
1 ~ 10 kb: 5 sec/ kb	1 ~ 10 kb: 5 sec/ kb
≥ 10 kb: 10 sec/ kb	≥ 10 kb: 10 sec/ kb

## Quality Control

Kit components are free of contaminating DNase and RNase. 2x repliQa HiFi ToughMix is functionally tested for amplification of a 4-kb fragment from a single-copy gene in a human genomic DNA.

Nuclease Assay:

DNase: DNase activity must be below the detectable limits of 100 pg DNase I equivalent as assayed using a fluorogenic substrate following a 1 hour incubation at 37°C with each kit component at 1x concentration.

RNase: RNase activity must be below the detectable limits of 1 pg RNase A equivalent as assayed using a fluorogenic substrate following a 1 hour incubation at 37°C with each kit component at 1x concentration.

4.1 kb PCR Functional Assay: Negative control must be free of visible product with a single band at ~4.1 kb visible from 35 cycles of PCR using 20 ng human genomic DNA.

## Limited Label Licenses

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this manual and for use with components contained in the kit only. Quantabio, LLC. grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the protocols provided with the product, this manual, and additional protocols available at [www.quantabio.com](http://www.quantabio.com). Some of these additional protocols have been provided by Quantabio product users. These protocols have not been thoroughly tested or optimized by Quantabio, LLC. Quantabio, LLC. neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, Quantabio, LLC. makes no warranty that this kit and/or its use(s) do not infringe the rights of third-parties.
3. This kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. Quantabio, LLC. specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. Quantabio, LLC. may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

©2021 Quantabio, LLC. 100 Cummings Center Suite 407J Beverly, MA 01915; Telephone number: 1-888-959-5165.

Quantabio products are manufactured in Beverly, Massachusetts, Frederick, Maryland and Hilden, Germany.

Intended for molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease.

## Trademarks

ToughMix is a registered trademark of Quantabio, LLC. TWEEN is a registered trademark of Croda International PLC

## Appendix A: Optimized PCR Product Cleanup

After PCR amplification, sample cleanup may be required to remove enzyme, dNTPs, and reaction buffer prior to sensitive downstream applications. The following reflects an optimized cleanup procedure using magnetic beads (i.e. sparQ PureMag Beads) to achieve high purity samples for further analysis.

When purifying highly concentrated or long DNA after amplification with repliQa HiFi ToughMix, abnormal pelleting of the beads may be observed, which can impede PCR product clean up. The following protocols can be used to improve bead clean up if this occurs.

### Additional reagents and materials that are not supplied

Proteinase K (concentration 20 mg/ml) for Option 1 ONLY

10% Tween® 20 for Option 2 ONLY

Purifications Beads

Purification beads are not included with this kit and must be purchased separately. This protocol has been validated using the sparQ PureMag Beads (Cat No. 95196-005, 95196-060, or 95196-450) from Quantabio.

### Reaction Cleanup with Proteinase K Digestion

1. Equilibrate the sparQ PureMag beads to RT.
2. If PCR reaction volume is less than 50  $\mu$ L, adjust the volume to 50  $\mu$ L using nuclease-free water or TE Buffer (10 mM Tris-HCl (pH 8.0), 0.1 mM EDTA).
3. Add 1  $\mu$ L Proteinase K to the sample, mix well and incubate at RT for at least 1 min.
4. Thoroughly vortex the bead slurry and add 51  $\mu$ L (1X) to the PCR reaction. Mix well by pipetting.
5. Incubate the mixture for 5 min at RT. Pellet the beads on a magnetic stand (e.g., DynaMag) and carefully discard the supernatant.
6. Keep the tube on the magnetic stand and wash the beads with 200  $\mu$ L of freshly prepared 80% ethanol.
7. Pellet the beads on the magnetic stand. Carefully discard the supernatant after the liquid is clear. Repeat the wash once.
8. Air-dry the beads on the magnetic stand for 5 min or until the beads are dry. Over-drying of beads may result in lower DNA recovery.
9. Resuspend the dried beads in 32.5  $\mu$ L or desired volume of 10 mM Tris-HCl, pH 8.0 and incubate for 2 min at RT. Pellet the beads on the magnetic stand. Carefully transfer 30  $\mu$ L of supernatant into a new tube. The sample can be stored at -20°C if not proceeding immediately to library quantification or other downstream processes.

### Reaction Cleanup with Tween® 20

1. Equilibrate the sparQ PureMag beads to RT.
2. If PCR reaction volume is less than 50 $\mu$ L, adjust the volume to 50  $\mu$ L using nuclease-free water or TE Buffer (10 mM Tris-HCl (pH 8.0), 0.1 mM EDTA).
3. Add 1  $\mu$ L 10% Tween® 20 to the sample and mix well.
4. Thoroughly vortex the beads slurry and add 51  $\mu$ L (1X) to the PCR reaction. Mix well by pipetting.
5. Incubate the mixture for 5 min at RT. Pellet the beads on a magnetic stand (e.g., DynaMag) and carefully discard the supernatant.
6. Keep the tube on the magnetic stand and wash the beads with 200  $\mu$ L of freshly prepared 80% ethanol.
7. Pellet the beads on the magnetic stand. Carefully discard the supernatant after the liquid is clear. Repeat the wash once.
8. Air-dry the beads on the magnetic stand for 5 min or until the beads are dry. Over-drying of beads may result in lower DNA recovery.
9. Resuspend the dried beads in 32.5  $\mu$ L or desired volume of 10 mM Tris-HCl, pH 8.0 and incubate for 2 min at RT. Pellet the beads on the magnetic stand. Carefully transfer 30  $\mu$ L of supernatant into a new tube. The sample can be stored at -20°C if not proceeding immediately to library quantification or other downstream processes.