

Single-tube solutions to streamline DNA library preparation and improve sequencing results while lowering costs

Yi Jin, Marissa Bolduc, David Bays, Shuhong Li, Eleanor Kolossovski, Brian Komorous, Hongbo Liu, and David Schuster
100 Cummings Center Suite 407J, Beverly, MA 01915


Introduction

The sparQ DNA Frag & Library Prep Kit optimizes the integration of enzymatic DNA fragmentation and library construction into a simple two-step workflow for Illumina® NGS platforms. An optimized mix of QuantaBio's engineered enzymes works in concert to combine tunable DNA fragmentation and polishing reactions minimizing over fragmentation and greatly simplifying the library prep process. Resulted 5'-phosphorylated and 3'-dA-tailed DNA fragments are suitable for direct ligation of sequencing adapters without an intervening cleanup step. The streamlined workflow can be completed in under 3 hours with minimal hands-on time accommodating DNA input amounts from 1 ng to 1 µg. The HiFi PCR Master Mix and Primer Mix allow for the unbiased amplification of DNA libraries. PCR-free workflow is enabled from 100 ng of DNA sample.

Features & Benefits

- ✓ Simple 2-step workflow employs a unique enzyme mix, safeguarding samples from over fragmentation
- ✓ Tunable and reproducible fragmentation profiles across a range of sample types
- ✓ Flexible generation of high quality libraries from 1 ng – 1 µg of input DNA
- ✓ PCR-free workflow enabled from 100 ng of input DNA
- ✓ Minimized bias across challenging regions for improved sequencing results

Library Prep Workflow




Workflow Comparisons


	NEBNext® Ultra™ II FS DNA Library Prep Kit	Nextera DNA Flex Library Prep Kit	KAPA HyperPlus Kit	sparQ DNA Frag & Library Prep Kit
Total / Hands-on Time	3 hr / 35 min	3 hr / 90 min*	2.5 hr / 35 min	2.5 hr / 30 min
Input	100 pg – 500 ng	1 – 500 ng	1 ng – 1 µg	1 ng – 1 µg**
Pros	Enzyme-based fragmentation	Tagmentation; Integrated normalization (100 – 500 ng)	Enzyme-based fragmentation	High library yield and sensitivity; Consistency; Flexibility
Cons	Low yield; Poor Reproducibility; GC Bias; Different workflows	Low yield; High duplication rate; Fixed fragment size; Different workflow for different inputs	Poor consistency and reproducibility	Optimization for different inputs
	<ul style="list-style-type: none"> Fragmentation, End Repair & A-tailing Ligation USER Excision Purification PCR Purification 	<ul style="list-style-type: none"> Tagmentation Purification PCR Purification 	<ul style="list-style-type: none"> Fragmentation End Repair & A-tailing Ligation Purification PCR Purification 	<ul style="list-style-type: none"> Fragmentation & End Polishing Ligation Purification PCR Purification

* vendor specified time ** down to 100 pg

Tunable Size Ranges for Varying DNA Inputs (1 ng – 1 µg)



Superior Efficiency & Yields




High Quality DNA Libraries


Fragmentation Method	1 ng Input DNA		100 ng Input DNA	
	Mapped Reads*	Duplication	Mapped Reads*	Duplication
sparQ DNA Frag & Lib Prep Kit	91.92%	0.07%	94.45%	0.04%
KAPA HyperPlus Kit	92.40%	0.08%	93.49%	0.03%
Covaris	92.96%	0.09%	93.61%	0.03%

sparQ DNA Frag & Library Prep Kit produces high quality libraries with high read alignment percentages and low duplication rates

Uniform Coverage Across Challenging Genomic Regions



sparQ Compared to Other Technologies



Conclusions

sparQ DNA Fragmentation & Library Prep Kit

- Enzyme Mix integrates fragmentation and DNA polishing reactions for a streamlined workflow
- Complete workflow in under 3 hours with minimal hands-on time
- Adjusting fragmentation time generates selective DNA sizes
- Low fragmentation and amplification bias enables more uniform genome coverage
- Excellent sequencing metrics
- Capable of PCR-free library prep workflow