

Date	14 May 2018
Issues	2

1 Important highlights from this release

Run reconstruction, rebuilding a run by uploading run data from the instrument, now works making connections to instruments even more dependable

Replicate statistics in the Cycling and Absolute Quantification analyses are now displayed in the reports.

The sample grouping settings in the Cycling and Absolute Quantification result tables are now saved to the run and project files. These settings are respected when run and project files are reopened.

2 All updates for this release

2.1 Improvement

- [CYCLERSW-1374](#) - Cycling and Absolute Quant Replicate statistics in report

2.2 Bug

- [CYCLERSW-1366](#) - Applications cannot reconstruct run