

| Product Information                       |                |  |  |
|---|----------------|--|--|
| PerfeCтa® qPCR FastMix®, UNG, Low<br>ROX™ |                |  |  |
| Part Number                               | 95078-05K      |  |  |
| Number of Reactions                       | 5000 Reactions |  |  |
| Reaction Size                             | 20 μL          |  |  |
| Storage Temperature                       | -25℃ to -15℃   |  |  |
| Lot Number                                | 66155747       |  |  |
| Reference Number                          | 073018         |  |  |
| Expiration Date                           | 07/31/2021     |  |  |

#### **Product Description:**

PerfeCTa qPCR FastMix, UNG is a 2X concentrated, ready-touse reaction cocktail that contains all components, except primers, probe(s), and template for real-time quantitative PCR systems that do not require an internal reference dye. The proprietary buffer and stabilizers have been specifically optimized to deliver maximum PCR efficiency, sensitivity, and robust fluorescent signal with TaqMan® or TaqMan MGB probe chemistry when using rapid PCR cycle times and reduced reaction volumes. This affords greater reagent economy and laboratory throughput on conventional or rapid ramp rate qPCR systems. The enhanced specificity of this FastMix suppresses cross-reactivity between homologous sequences, improving detection and discrimination in SNP applications. A key component of this FastMix is AccuFast™ Tag DNA polymerase. This hot-start Tag contains a proprietary mixture of monoclonal antibodies that bind to the polymerase and keep it inactive prior to the initial PCR denaturation step (> 48 hours at room temperature). Similar to our AccuStart™ Taq DNA polymerase, these antibodies are irreversibly inactivated

during the initial PCR denaturation step. However, unlike other antibody hot-start polymerases, activation of AccuFast Taq is instantaneous at 95°C. Rapid recovery of fully active, unmodified Taq DNA polymerase is critical for efficient extension kinetics. Replication of fragments up to 200 bp is complete in less than 20s at 60°C. Additionally, the dNTP mix in this Fast Mix contains dUTP in place of dTTP. Inclusion of uracil-N-glycosylase (UNG) prevents amplification of carry-over contamination from previous dU-containing PCRs.

## **Component Part Numbers:**

84082 PerfeCTa gPCR FastMix, UNG, Low ROX 50 mL

| Product Specifications |                  |       |       |
|------------------------|------------------|-------|-------|
| 95078                  |                  |       |       |
| Assay                  | Functional Assay | DNase | RNase |
| Result                 | Pass             | Pass  | Pass  |

### **Quality Control Analysis and Specifications:**

#### **Nuclease Assay:**

**DNase:** DNase activity must be below the detectable limits of 100 pg DNase I equivalent as assayed using a fluorogenic substrate following a 1 hour incubation at 37°C with each kit component at 1X concentration.

**RNase:** RNase activity must be below the detectable limits of 1 pg RNase A equivalent as assayed using a fluorogenic substrate following a 1 hour incubation at 37°C with each kit component at 1X concentration.

**qPCR B-actin Plasmid DNA Functional Assay:** Detection of B-actin from 10 copies to 1 x  $10^7$  copies. The Cq standard curve analysis must have a coefficient of determination ( $r^2$ )  $\geq 0.990$  with a slope between -3.20 to -3.70

# Limitations of Use

Quantabio and Ultraplex are registered trademarks of QIAGEN Beverly, Inc. Quanta Biosciences, qScript, Geltrack, Tough Mix, PerfeCta, and Fastmix are registered trademarks of Quanta BioSciences Inc. Extracta, AccuStart, AccuMelt, and Accuvue are trademarks of Quanta BioSciences Inc. Applied Biosystems, StepOne, StepOnePlus and ROX are trademarks of Thermo Fisher Scientific and or its subsidiaries. Please contact QIAGEN-Beverly for more information.

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for a dministration to humans or a nimals. SDS sheets relevant to this product are available upon request.