

Product Information	
repliQa™ HiFi Assembly Mix	
Part Number	95190-050
Number of Reactions	50 Reactions*
Reaction Size	20 µL
Storage Temperature	-25°C to -15°C
Lot Number	025928
Reference Number	071717, 021518, D5061818
Expiration Date	02/28/2019

\*Variable depending on individual usage of HiFi Assembly Mix

**Product Description:**

The repliQa HiFi Assembly Mix simplifies the construction of recombinant DNA through the simultaneous and seamless assembly of multiple DNA fragments possessing terminal regions of sequence overlap in a single, isothermal reaction. Similar in principle to the Gibson Assembly® Method, the high efficiency repliQa HiFi Assembly Mix is ideal for a range of genetic engineering applications including routine molecular cloning, site-directed mutagenesis, assembly of large constructs for synthetic biology applications, and the construction of diverse sequence libraries for directed evolution studies. The concentrated, two-component format allows flexibility in design of assembly reactions and compatibility with less concentrated DNA samples. The repliQa Mix has been optimized for use with a total input quantity of DNA fragments in the range of 0.03 to 0.5 pmols. The assembly of up to six DNA fragments is recommended, though the repliQa Mix has been successfully used for more complex assemblies.

**Component Part Numbers:**

**84616 repliQa HiFi Assembly Enzyme Mix 0.10ml**

**84617 repliQa 10X Assembly Reaction Buffer 0.50ml**

**84619 repliQa DpnI Restriction Enzyme 0.25ml**

Product Specifications	
95190	
Assay	DNA Assembly Quantification Assay
Specification	≥15%

**Quality Control Analysis and Specifications:**

The repliQa HiFi Assembly Mix is functionally tested for assembly of three 1-kb PCR fragments and quantitation of the 3 kb product.

Enzyme components were tested prior to assembly and found free of contaminating endonucleases and exonucleases. Enzyme purity was >95% as determined by SDS-PAGE and negligible *E.coli* genomic DNA contamination was confirmed by qPCR.

**Limitations of Use**

QuantaBio and Ultrplex are registered trademarks of QIAGEN Beverly, Inc. Quanta Biosciences, qScript, Geltrack, ToughMix, PerfeCta, and Fastmix are registered trademarks of Quanta BioSciences Inc. Extracta, AccuStart, AccuMelt, and Accuvue are trademarks of Quanta BioSciences Inc. Applied Biosystems, StepOne, StepOnePlus and ROX are trademarks of Thermo Fisher Scientific and or its subsidiaries. Please contact QIAGEN-Beverly for more information.

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. SDS sheets relevant to this product are available upon request.

100 Cummings Center, Suite 407J, Beverly, MA 01915 • Ph (888) 927-7027 • Fax (978) 867-5724 • [www.QuantaBio.com](http://www.QuantaBio.com) FMWI016.2 Rev 01