# Quantabio

Product Information			
PerfeCTa SYBR <sup>®</sup> Green SuperMix low ROX <sup>®</sup>			
95056-500			
500 reactions			
50µL			
-25°C to -15°C			
66198623			
012122			
12/31/2022			

## **Product Description:**

PerfeCTa SYBR Green SuperMix low Rox is a user-friendly, 2X concentrated reaction mix that simplifies setup and reduces errors with optimized reference dye and pre-blended AccuVue plate loading dye for visual confirmation of reagent addition and mixing. This proprietary buffer technology stabilizes a high concentration of SYBR Green I dye to ensure maximum optical signal with low abundance or small targets (such as microRNA). Successful detection with a non-specific, dsDNA intercalating dye requires precise target amplification as off-target primer elongation will contribute to overall fluorescent signal and lead to over-reported relative abundance values. This reagent is powered by a highlyprocessive, ultra-pure Taq DNA polymerase mutant with stringent, ultra-pure AccuStart™II antibody hot start technology that allows ambient room-temperature setup and maximal enzyme kinetics after rapid, irreversible denaturation at 95°C.

## Component Part Numbers: 84020 PerfeCTa SYBR Green SuperMix Low ROX 1.25mL

Product Specifications			
95056			
Assay	qPCR ß-Actin Plasmid DNA Functional Assay	DNase	RNase
Result	Pass	Pass	Pass

## **Quality Control Analysis and Specifications:**

## Nuclease Assay:

**DNase:** DNase activity must be below the detectable limits of 100 pg DNase I equivalent as assayed using a fluorogenic substrate following a 1 hour incubation at 37°C with each kit component at 1X concentration.

**RNase:** RNase activity must be below the detectable limits of 1 pg RNase A equivalent as assayed using a fluorogenic substrate following a 1 hour incubation at 37°C with each kit component at 1X concentration.

**qPCR β-Actin Plasmid DNA Functional Assay for SYBR Green SuperMix and SYBR Green FastMix:** Detection of β-actin from 10 copies to  $1 \times 10^7$  copies. The Coefficient of Determination ( $r^2$ )  $\ge 0.990$  with a slope analysis between -3.20 and -3.65

### Limitations of Use

Quantabio and Ultraplex are registered trademarks of QIAGEN Beverly, Inc. Quanta Biosciences, qScript, Geltrack, ToughMix, PerfeCta, and Fastmix are registered trademarks of Quanta BioSciences Inc. Extracta, AccuStart, AccuMelt, and Accuvue are trademarks of Quanta BioSciences Inc. Applied Biosystems, StepOne, StepOnePlus and ROX are trademarks of Thermo Fisher Scientific and or its subsidiaries. Please contact QIAGEN-Beverly for more information.

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. SDS sheets relevant to this product are available upon request. 100 Cummings Center, Suite 407J, Beverly, MA 01915 • Ph (888) 927-7027 • Fax (978) 867-5724 • <u>www.QuantaBio.com</u> FMWI016.2 Rev 01