

| Product Information | |
|-----------------------|----------------|
| repliQa HiFi ToughMix | |
| Part Number | 95200-025 |
| Number of Reactions | 25 reactions |
| Reaction Size | 50 uL |
| Storage Temperature | -25°C to -15°C |
| Lot Number | 66228309 |
| Reference Number | 153300 |
| Expiration Date | 10/31/2024 |

Product Description:

The repliQa HiFi ToughMix is a 2x, ready-to-use solution that contains all the components for high fidelity PCR amplification, including a genetically modified DNA polymerase coupled with hot start antibodies. This unique, next generation master mix provides 90x higher fidelity compared to Taq, while reducing time to PCR results by 2-3x. The extreme speed is enabled by extension times as fast as 1-10 kb/sec depending on target length. The enzyme is coupled with the industry leading ToughMix which is tolerant to a wide variety of inhibitors making it suitable for routine PCR, cloning, amplicon sequencing and site directed mutagenesis.

Component Part Numbers:

84640 repliQa HiFi ToughMix 0.625 mL

| Product Specifications | |
|------------------------|--|
| 95200 | |
| Assay | repliQa HiFi ToughMix Functional Assay |
| Result | Pass |

Quality Control Analysis and Specifications:

The functionality of the repliQa HiFi ToughMix is evaluated in a PCR reaction to produce a single 4.1kb product from 20 ng of Human Genomic DNA, while also ensuring the negative control is free of visible product.

Limitations of Use

QuantaBio and Ultrplex are registered trademarks of QIAGEN Beverly, Inc. Quanta Biosciences, qScript, Geltrack, ToughMix, PerfeCta, and Fastmix are registered trademarks of Quanta BioSciences Inc. Extracta, AccuStart, AccuMelt, and Accuvue are trademarks of Quanta BioSciences Inc. Applied Biosystems, StepOne, StepOnePlus and ROX are trademarks of Thermo Fisher Scientific and or its subsidiaries. Please contact QIAGEN-Beverly for more information.

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. SDS sheets relevant to this product are available upon request.

100 Cummings Center, Suite 407J, Beverly, MA 01915 • Ph (888) 927-7027 • Fax (978) 867-5724 • www.enzymatics.com FMWI016.2 Rev 01