

Product Information	
sparQ UDI Adapters (1-24)	
Part Number	95211-024
Number of Reactions	24 reactions
Reaction Size	N/A
Storage Temperature	-25 to -15°C
Lot Number	66296522
Reference Number	178036228
Expiration Date	N/A

Product Description:

sparQ UDI Adapters (1-24) are unique dual-indexed adapters designed to be used on Illumina sequencing platforms. The use of UDI's minimizes index hopping events while also enabling greater multiplexing on flow cells.

Well Position	Barcode Number	Barcode Sequence
1	QDIY-001	TGAACGTTGTATGGCCGACT
2	QDIY-002	ACCAGACTTGCGATGAGCAC
3	QDIY-003	ACTGGCGAACGATAAGTCGA
4	QDIY-004	GCGTTAGGCATCACGCCTTG
5	QDIY-005	TTATCGGCCTAGGAACACAA
6	QDIY-006	GAGGTATAAGCTCAGTAGGC
7	QDIY-007	TCAAGGATTCTGAAGTGCCTG
8	QDIY-008	CGAACCGAGATCTCTCGCCT
9	QDIY-009	GAGCCAAGTTAGGCACCTTC
10	QDIY-010	AAGGCCGTAGCTGTTGGTAA
11	QDIY-011	TTAGAGAAGCGCTGGTACCT
12	QDIY-012	TCTAAGACCATAAGGAGCGG

Well Position	Barcode Number	Barcode Sequence
13	QDIY-013	TGTAACCACTAATCGCTCCA
14	QDIY-014	CCGACACAAGCTCCTAATTG
15	QDIY-015	CTCTGATGGCGCCTCATAAT
16	QDIY-016	CGGCCTGTTATGTATTGAGC
17	QDIY-017	TGCATAGCTTAGCCATAACA
18	QDIY-018	AACCTTCTCGCCACAAGTGG
19	QDIY-019	AAGAGATCACGTTATCACAC
20	QDIY-020	GCCTGAAGGATACCGTTCTT
21	QDIY-021	ATTGTGCCTTAGGCGTTAGG
22	QDIY-022	TCCTCTACCGCCGTAACGTC
23	QDIY-023	TACCATGAACGTAATAGCCA
24	QDIY-024	CATTGGCAGATAGCGCCGAT

Limitations of Use

Quantabio, UltraPlex, qScript, GelTrack, ToughMix, PerfeCta, and FastMix are registered trademarks of Quantabio, LLC. Applied Biosystems, StepOne, StepOnePlus and ROX are trademarks of Thermo Fisher Scientific and or its subsidiaries. Please contact Quantabio for more information.

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. SDS sheets relevant to this product are available upon request.

100 Cummings Center, Suite 407J, Beverly, MA 01915 • Phone +1 888-959-5165 • www.quantabio.com FMWI-016.2 Rev.03