

Product Information	
<b>sparQ HiFi PCR Master Mix 50R</b>	
Part Number	95192-050
Number of Reactions	50 Reactions
Reaction Size	50 µL
Storage Temperature	-25°C to -15°C
Lot Number	66175721
Reference Number	011821, 275754821, 9556003
Expiration Date	10/31/2022

**Product Description:**

A high efficiency, high-fidelity, and low bias 2X PCR master mix for NGS library amplification.

**Component Part Numbers:**

84523 HiFi PCR Master Mix (2X), 1.25mL  
 84526 Primer Mix, 0.15mL  
 84527 HiFi Enhancer, 0.150

Product Specifications	
95192	
Assay	HiFi PCR Master Mix Functional Assay
Specification	Functional

**Quality Control Analysis and Specifications:**

**HiFi PCR Master Mix Functional Assay:** Quality of the 2X HiFi PCR Master Mix is tested functionally by amplification of a DNA library prepared from mixed bacterial genomic DNA with GC-content of 10-80%. The differences in library yield and profile among different lots must not exceed 15%. Sequencing of the amplified library must yield mapped reads >90% and normalized coverage between 0.7 and 1.3 across the full GC spectrum.

**Limitations of Use**

QuantaBio and Ultrplex are registered trademarks of QIAGEN Beverly, Inc. Quanta Biosciences, qScript, Geltrack, ToughMix, PerfeCta, and Fastmix are registered trademarks of Quanta BioSciences Inc. Extracta, AccuStart, AccuMelt, and Accuvue are trademarks of Quanta BioSciences Inc. Applied Biosystems, StepOne, StepOnePlus and ROX are trademarks of Thermo Fisher Scientific and or its subsidiaries. Please contact QIAGEN-Beverly for more information.

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. SDS sheets relevant to this product are available upon request.