

Product Information	
<b>Twist Library Preparation Kit 1</b>	
QIAGEN Part Number	CMK0014
Twist P/N	100572
Unit Size	96 reactions
Storage Temperature	-25°C to -15°C
Twist Lot Number	3200000197
Twist Expiry Date	05/2022
Kit Manufacturing Date	12/14/2020
Reference Number	66172621

**Product Description:**

Twist Library Preparation Kit 1

Kit Components						
Description	QIAGEN P/N	QIAGEN Lot Number	QIAGEN Expiry	Twist P/N	Twist Lot Number	Twist Expiry
5x Fragmentation Enzyme	Y9410C-4	111820R	05/2022	100579	3200000196	05/2022
10x Fragmentation Buffer	B0330C-4	111220	05/2022	100580	3200000195	05/2022
DNA Ligation Mix	L6030C-11	84120920	11/2022	100581	3200000198	11/2022
DNA Ligation Buffer	B9020C-6	082020	08/2023	100582	3200000194	08/2023
Amplification Primers, ILMN	CM0171C-1	N/A	N/A	100583	3200000192	12/2022

Product Specifications						
L6030						
Assay	SDS Purity	Specific Activity	SS Exonuclease	DS Exonuclease	DS Endonuclease	E. coli DNA Contamination
Units Tested	n/a	n/a	6,000 U	6,000 U	6,000 U	6,000 U
Specification	>99%	300,000 U/mg	< 1.0 % Released	< 1.0 % Released	No Conversion	< 10 copies

**Quality Control Analysis:**

Enzyme components were tested prior to assembly and found free of contaminating endonucleases and exonucleases. Enzyme purity was >95% as determined by SDS-PAGE and negligible *E.coli* genomic DNA contamination was confirmed by qPCR. Specific activity was verified for each enzyme.

**L6030:**

**Unit Activity** is measured using a 2-fold serial dilution method. Dilutions of enzyme batch were made in 1X DNA Ligase Reaction Buffer and added to 20 µL reactions containing double stranded DNA fragments and 1X DNA Ligase Reaction Buffer. Reactions are incubated for 30 minutes at 23°C, stopped, and analyzed on a 1% agarose gel stained with ethidium bromide.

**Protein Concentration** is determined by OD<sub>280</sub> absorbance.

**Limitations of Use**

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. MSDS sheets relevant to this product are available upon request.

**Physical Purity** is evaluated by SDS-PAGE of concentrated and diluted enzyme solutions followed by silver stain detection. Purity is assessed by comparing the aggregate mass of contaminant bands in the concentrated sample to the mass of the protein of interest band in the diluted sample.

**Single-Stranded Exonuclease** is determined in a 50 µL reaction containing a radiolabeled single-stranded DNA substrate and 10 µL of enzyme solution incubated for 4 hours at 37°C.

**Double-Stranded Exonuclease** is determined in a 50 µL reaction containing a radiolabeled double-stranded DNA substrate and 10 µL of enzyme solution incubated for 4 hours at 37°C.

**Double-Stranded Endonuclease** is determined in a 50 µL reaction containing 0.5 µg of plasmid DNA and 10 µL of enzyme solution incubated for 4 hours at 37°C.

**E.coli 16S rDNA Contamination** is evaluated using 5 µL replicate samples of enzyme solution denatured and screened in a TaqMan qPCR assay for the presence of contaminating *E.coli* genomic DNA using oligonucleotide primers corresponding to the 16S rRNA locus.

**Y9410:**

**The functionality of the 5X WGS Fragmentation Mix** is evaluated by performing library construction, PCR/qPCR and sequencing.

Analysis	Acceptance Criteria
Fragmentation of 100ng human genomic DNA	<ul style="list-style-type: none"> <li>300 bp +/-30 bp peak size of DNA fragment profile is expected at tested condition</li> <li>Average fragment size from test lots must be within 15% of results from the reference lot</li> </ul>
Quantification of library	Final library concentration must be >60nM
Sequencing	<ul style="list-style-type: none"> <li>Normalized coverage between 0.7-1.3 for tested GC-contents (10– 80%)</li> <li>mapped reads must be &gt;90%</li> <li>Q score of the sequencing run &gt;85%</li> </ul>

**Limitations of Use**

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. MSDS sheets relevant to this product are available upon request.