

Product Information	
<b>Twist Library Preparation Kit 1: Mechanical Fragmentation, 96 Samples</b>	
QIAGEN Part Number	CMK0022-96
Twist P/N	100876
Unit Size	96 reactions
Storage Temperature	-25°C to -15°C
Twist Lot Number	3200000283
Twist Expiry Date	04/01/2023
Kit Manufacturing Date	10/14/2021
Reference Number	66191649

**Product Description:**

Twist Library Preparation Kit 1: Mechanical Fragmentation

Kit Components						
Description	QIAGEN P/N	QIAGEN Lot Number	QIAGEN Expiry	Twist P/N	Twist Lot Number	Twist Expiry
5X ER/A-Tailing Enzyme Mix	Y9420C-6	083021	08/2023	101222	3200000278	08/2023
10X ERA Buffer	B9420C-6	051421	05/2023	101223	3200000279	05/2023
DNA Ligation Mix	L6030C-11	110092021	09/2023	100581	3200000280	09/2023
DNA Ligation Buffer	B9020C-6	081721	08/2024	100582	3200000284	08/2024
Amplification Primers, ILMN	CM0171C-1	N/A	N/A	100583	3200000282	04/2023

Product Specifications	
Y9240	
Assay	ER/A Tailing Enzyme Mix Functional Assay
Specification	Functional

Product Specifications						
L6030						
Assay	SDS Purity	Specific Activity	SS Exonuclease	DS Exonuclease	DS Endonuclease	E. coli DNA Contamination
Units Tested	n/a	n/a	6,000 U	6,000 U	6,000 U	6,000 U
Specification	>99%	300,000 U/mg	< 1.0 % Released	< 1.0 % Released	No Conversion	< 10 copies

**Quality Control Analysis:**

Enzyme components were tested prior to assembly and found free of contaminating endonucleases and exonucleases. Enzyme purity was >95% as determined by SDS-PAGE and negligible *E. coli* genomic DNA contamination was confirmed by qPCR. Specific activity was verified for each enzyme.

**Limitations of Use**

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. MSDS sheets relevant to this product are available upon request.

**L6030:**

**Unit Activity** is measured using a 2-fold serial dilution method. Dilutions of enzyme batch were made in 1X DNA Ligase Reaction Buffer and added to 20 µL reactions containing double stranded DNA fragments and 1X DNA Ligase Reaction Buffer. Reactions are incubated for 30 minutes at 23°C, stopped, and analyzed on a 1% agarose gel stained with ethidium bromide.

**Protein Concentration** is determined by OD<sub>280</sub> absorbance.

**Physical Purity** is evaluated by SDS-PAGE of concentrated and diluted enzyme solutions followed by silver stain detection. Purity is assessed by comparing the aggregate mass of contaminant bands in the concentrated sample to the mass of the protein of interest band in the diluted sample.

**Single-Stranded Exonuclease** is determined in a 50 µL reaction containing a radiolabeled single-stranded DNA substrate and 10 µL of enzyme solution incubated for 4 hours at 37°C.

**Double-Stranded Exonuclease** is determined in a 50 µL reaction containing a radiolabeled double-stranded DNA substrate and 10 µL of enzyme solution incubated for 4 hours at 37°C.

**Double-Stranded Endonuclease** is determined in a 50 µL reaction containing 0.5 µg of plasmid DNA and 10 µL of enzyme solution incubated for 4 hours at 37°C.

**E.coli 16S rDNA Contamination** is evaluated using 5 µL replicate samples of enzyme solution denatured and screened in a TaqMan qPCR assay for the presence of contaminating *E.coli* genomic DNA using oligonucleotide primers corresponding to the 16S rRNA locus.

**Y9420:**

**ER/A Tailing Enzyme Mix Functional Assay:** QC Library length must be within 15% of the reference library length. Concentration of the QC library generated from 100 ng input DNA (average ~300 bp fragments) is >60 nM with mapped reads > 90%. For QC library, normalized coverage should be within 0.7 to 1.3 for most of the genome (10% - 80% GC content).

---

**Limitations of Use**

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. MSDS sheets relevant to this product are available upon request.