

Product Information	
<b>5PRIME HotMaster® Taq DNA Polymerase</b>	
Part Number	2200300
Number of Units	100 Units
Storage Temperature	-25°C to -15°C
Lot Number	022434
Expiration Date	01/31/2019

**Product Description:**

Innovative Hot-Start/Cold-Stop technology for highly specific Hot-Start PCR without enzyme activation

**Features and Benefits**

- Highly specific amplification
- Minimal optimization with self-adjusting Mg<sup>2+</sup> buffer technology
- No enzyme activation - excellent for fast PCR protocols
- Convenient reaction set-up at room temperature

**Component Part Numbers:**

84508 HotMaster Taq DNA Polymerase 0.02 mL

84509 10X HotMaster Taq Buffer 1.8mL

Product Specifications	
<b>84508 HotMaster Taq DNA Polymerase</b>	
Assay	1.1 Kb PCR Fragment Assay
Result	Pass

Product Specifications	
<b>84509 10X HotMaster Taq Buffer</b>	
Assay	Conductivity
Result	39.48 – 45.32 mS/cm

**Quality Control Analysis and Specifications:**

**1.1 Kb PCR Fragment Assay:** The major PCR product should be a 1.1 Kb PCR Fragment

**Limitations of Use**

HotMaster is a registered trademarks of QIAGEN Beverly, Inc Please contact QIAGEN-Beverly for more information.

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. SDS sheets relevant to this product are available upon request.

100 Cummings Center, Suite 407J, Beverly, MA 01915 • Ph (888) 927-7027 • Fax (978) 867-5724 • [www.5prime.com](http://www.5prime.com) FMW1016.2 Rev 01