

Product Information	
<b>PerfeCTa® SYBR® Green SuperMix</b>	
Part Number	95054-500
Number of Reactions	500 reactions
Reaction Size	50 µL
Storage Temperature	-25°C to -15°C
Lot Number	021931
Reference Number	121916
Expiration Date	12/31/2019

**Product Description:**

PerfeCTa SYBR Green SuperMix is a user-friendly, 2X concentrated reaction mix that simplifies setup and reduces errors with optimized reference dye and pre-blended AccuVue plate loading dye for visual confirmation of reagent addition and mixing. This proprietary buffer technology stabilizes a high concentration of SYBR Green I dye to ensure maximum optical signal with low abundance or small targets (such as microRNA). Successful detection with a non-specific, dsDNA intercalating dye requires precise target amplification as off-target primer elongation will contribute to overall fluorescent signal and lead to over-reported relative abundance values. This reagent is powered by a highly-processive, ultra-pure Taq DNA polymerase mutant with stringent, ultra-pure AccuStart™II antibody hot start technology that allows ambient room-temperature setup and maximal enzyme kinetics after rapid, irreversible denaturation at 95°C.

**Component Part Numbers:**

84016 PerfeCTa SYBR Green SuperMix, 1.25mL

Product Specifications			
95054			
Assay	qPCR β actin Plasmid DNA Functional Assay for SYBR Green SuperMix	DNase	RNase
Result	Pass	Pass	Pass

**Quality Control Analysis and Specifications:**

**Nuclease Assay:**

**DNase:** DNase activity must be below the detectable limits of 100 pg DNase I equivalent as assayed using a fluorogenic substrate following a 1 hour incubation at 37°C with each kit component at 1X concentration.

**RNase:** RNase activity must be below the detectable limits of 1 pg RNase A equivalent as assayed using a fluorogenic substrate following a 1 hour incubation at 37°C with each kit component at 1X concentration.

**qPCR β actin Plasmid DNA Functional Assay for SYBR Green SuperMix:** Real-time PCR detection of log-fold serial dilutions of a control DNA from 10 copies to 1 x 10<sup>7</sup> copies. Cq standard curve analysis must have coefficient of determination (r<sup>2</sup>) ≥0.990 with a slope between -3.20 and -3.65.

**Limitations of Use**

QuantaBio and Ultrplex are registered trademarks of QIAGEN Beverly, Inc. Quanta Biosciences, qScript, Geltrack, ToughMix, PerfeCTa, and Fastmix are registered trademarks of Quanta BioSciences Inc. Extracta, AccuStart, AccuMelt, and Accuvue are trademarks of Quanta BioSciences Inc. Applied Biosystems, StepOne, StepOnePlus and ROX are trademarks of Thermo Fisher Scientific and or its subsidiaries. Please contact QIAGEN-Beverly for more information.

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. SDS sheets relevant to this product are available upon request.