

Product Information	
qScript® Flex cDNA Kit	
Part Number	95049-025
Number of Reactions	25 Reactions
Reaction Size	20 µL
Storage Temperature	-25°C to -15°C
Lot Number	66203737
Reference Number	101921, 062121, 021721, 081821, 071521
Expiration Date	06/30/2024

Product Description:

The qScript Flex cDNA Synthesis Kit is configured in separate components for supporting multiple RT priming strategies. The kit provides optimized, flexible components for priming with oligo-dT(20), random hexamers, gene-specific primer (GSP), or any combination thereof. This unique 5X master mix of buffer, magnesium, stabilizers and dNTPs simplifies reaction assembly and ensures robust and reproducible synthesis of first strand product from 10 pg to 4µg of total RNA or purified polyA+ RNA template. The resulting cDNA product is directly compatible with current real-time RT-PCR methods or end-point RT-PCR. Length of cDNA product is dependent upon priming strategy and quality of the RNA template. Oligo dT or GSP can be used for two-step RT-PCR of RNA targets up to 12 kb. Random primer is suitable for RT-PCR of RNA targets less than 1 kb. Any of the three priming methods are suitable for real-time quantitative RT PCR (RT-qPCR).

Component Part Numbers:

- 84001 qScript RT 0.025 mL
- 84039 5X qScript Flex Reaction Mix 0.10 mL
- 84040 Oligo dT 0.050 mL
- 84041 Random Primer 0.050 mL
- 84042 GSP Enhancer 0.050 mL
- 84007 Nuclease Free Water 1.50 mL

Product Specifications				
95049				
Assay	Functional PCR Assay for Flex Kit	β-actin SYBR Green qRT-PCR Assay for qScript Reverse Transcriptase	DNase	RNase
Result	Pass	Pass	Pass	Pass

Quality Control Analysis and Specifications:

Functional PCR Assay for Flex Kit Reaction Mix (5X): Detection of β-actin mRNA from 100 ng to 100 fg of total RNA. Coefficient of determination (R^2) ≥ 0.990 with a slope analysis between -3.20 and -3.70

Single bands visible at 6kb and 8kb from 35 cycles of PCR using 100 ng Oligo-dT cDNA

β-actin SYBR Green qRT-PCR Assay for qScript Reverse Transcriptase: Detection of β-actin mRNA from 100ng to 100fg of total RNA. Coefficient of determination (R^2) ≥ 0.990 with a slope analysis between -3.20 and -3.70

Nuclease Assay:

DNase: DNase activity must be below the detectable limits of 100 pg DNase I equivalent as assayed using a fluorogenic substrate following a 1 hour incubation at 37°C with each kit component at 1X concentration.

RNase: RNase activity must be below the detectable limits of 1 pg RNase A equivalent as assayed using a fluorogenic substrate following a 1 hour incubation at 37°C with each kit component at 1X concentration.

Limitations of Use

QuantaBio and Ultralex are registered trademarks of QIAGEN Beverly, Inc. Quanta Biosciences, qScript, Geltrack, ToughMix, PerfeCta, and Fastmix are registered trademarks of Quanta BioSciences Inc. Extracta, AccuStart, AccuMelt, and Accuvue are trademarks of Quanta BioSciences Inc. Applied Biosystems, StepOne, StepOnePlus and ROX are trademarks of Thermo Fisher Scientific and or its subsidiaries. Please contact QIAGEN-Beverly for more information.

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. SDS sheets relevant to this product are available upon request.



Product Specifications (Product Number) Rev (Number)

Limitations of Use

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. MSDS sheets relevant to this product are available upon request.