

Product Information	
Extracta™ DBS	
Part Number	95171-010
Volume	10.0 mL
Storage Temperature	2°C to 8°C
Lot Number	66240065
Reference Number	081723
Expiration Date	08/31/2025

**Product Description:**

Extracta DBS is a ready-to-use DNA extraction reagent for rapid and efficient recovery of PCR-ready DNA from dried blood spots (DBS) on Guthrie cards or Whatman 903 filter paper. This patented single-solution process produces DNA that is suitable for a variety of downstream applications including real-time qPCR or template generation for Next Generation Sequencing (NGS) or Sanger sequencing (1-3). DNA eluates are substantially free of PCR inhibitors and compatible with a variety of end-point PCR or real-time qPCR reagents. Application of Extracta DBS with PerfeCta® qPCR ToughMix® or PerfeCta® MultiPlex qPCR ToughMix enables accurate and reproducible quantification of DNA sequences in blood using TaqMan® hydrolysis probe real-time qPCR.

Product Specifications	
95171	
Assay	pH
Result	11.80 ± 0.2 @22°C

**Limitations of Use**

QuantaBio, UltraPlex, qScript, GelTrack, ToughMix, PerfeCta, and FastMix are registered trademarks of QuantaBio, LLC. Applied Biosystems, StepOne, StepOnePlus and ROX are trademarks of Thermo Fisher Scientific and or its subsidiaries. Please contact QuantaBio for more information.

This product was developed, manufactured, and sold for *in vitro* use only. The product is not suitable for administration to humans or animals. SDS sheets relevant to this product are available upon request.

100 Cummings Center, Suite 407J, Beverly, MA 01915 • Phone +1 888-959-5165 • [www.quantabio.com](http://www.quantabio.com) FMWI-016.2 Rev.03